

ALLURA[®]

Vaginal Stent



Intended to maintain the vaginal canal following radiological treatment or surgical procedures, ALLURA[®] assists in restoring, enlarging or creating a vagina.

- Made of an implant grade silicone shell and foam core, which is designed to expand and maintain the shape of the vagina without causing the patient additional discomfort
- Incorporates a locking ball check valve. This is able to be connected directly to a filling syringe, allowing the stent to be filled with sterile saline or air
- Integrated drainage tube is capable of being connected directly to a syringe to facilitate surgical site irrigation

Allura[®] is a single patient-use device, supplied sterile and intended for up to 90 days of use.

Allura[®] is offered in seven different standard sizes. PMT[®] can also create custom sized options based on your patient's unique requirements. To learn about custom PMT[®] products, reach out to your local PMT[®] representative.

MODEL NUMBER	Diam. X Length (CM)	MAX FILL (CC)
3321-09	3 X 9.5	64
3321-10	3.5 X 10	91
3321-11	4 X 11	132
3321-12	4 X 12	145
3321-14	4 X 14	170
3321-15	4 X 15	182
3321-16	4 X 16.5	200

Suture tabs are an available option for this product. When ordering, add "-T" to the end of the model number.

PMT[®]
PMT CORPORATION
QUALITY, INNOVATION AND CUSTOMER DRIVEN

PATENT PENDING
PMT, ALLURA and MANKIND are registered trademarks of PMT Corporation.
© Copyright 2020 PMT Corporation. All Rights Reserved.
PN 22912 Rev. 0 09/2020

Call Toll Free, 24 Hours a Day, 365 Days a Year
Phone: 1-800-MANKIND | Fax: 952-470-0865 | www.pmtcorp.com
info@pmtcorp.com