

# PMT<sup>®</sup> INTEGRA<sup>®</sup>

B R E A S T   E X P A N D E R S

## New Dimensions in Breast Expansion



Integra<sup>®</sup> breast expanders provide unmatched design features and characteristics for optimal results. With the two stage expansion models the Integra breast expanders provide greater lower pole projection to create a more anatomical expansion. Combined with the enhanced surface and optional suture tabs the Integra breast expander will provide you and your patient with the results desired. The Integra breast expanders are available in over 50 standard shapes and sizes but if your patient has unique requirements we are able to design a specialty breast expander to fill your needs. Please contact your PMT product specialist for more information.

# PRODUCT INFORMATION


# INTEGRA®

B R E A S T E X P A N D E R S



## PMT Breast Expanders Provide:

- » Lower pole projection that expands for a more natural shape.
- » A port which is easily identifiable through the skin.
- » A unique Enhanced™ surface which provides superior capillary adhesion.
- » A low profile design for easy placement and positioning.



	SIZE (CM) A x B	FILL (CC)	PROJECTION (CM)	SMOOTH #	ENHANCED #
Teardrop	10 x 13	500	8.5	3612-22-2i	3612es-22-2i
	11 x 14	600	10	3612-23-2i	3612es-23-2i
	12 x 16	850	11	3612-24-2i	3612es-24-2i
	14 x 18	1200	13	3612-26-2i	3612es-26-2i
Rounded Teardrop	10 x 11	400	8.5	3612-60-2i	3612es-60-2i
	11 x 12	550	9.5	3612-61-2i	3612es-61-2i
	11.5 x 12	650	10	3612-62-2i	3612es-62-2i
	12.5 x 13.5	850	11	3612-63-2i	3612es-63-2i
	14 x 14	1200	12	3612-64-2i	3612es-64-2i
Wide Base Round	11 x 11	500	9.5	3611-38-2i	3611es-38-2i
	13 x 13	900	12.5	3611-27-2i	3611es-27-2i
Wide Base Oval	10 x 12	450	9	3611-204-2i	3611es-204-2i
	11 x 13	650	10	3611-200-2i	3611es-200-2i
	12 x 14	850	11	3611-552-2i	3611es-552-2i
	13 x 15	1000	11	3611-199-2i	3611es-199-2i
	14 x 16	1200	12.5	3611-553-2i	3611es-553-2i

\* All models are available with Internal or Remote ports and Suture tabs. Additional models available, contact a PMT sales representative.

▶ **Local sales representatives are available 24 hours a day, 7 days a week, 365 days a year to assist you. Contact PMT® Corporation for more information.**

FDA REGISTERED AND APPROVED 

These products are protected by one or more of the following U.S. patents: 5,525,275; 4,685,447; 5,964,803; 5,456,716

PMT, ACCUSPAN, INTEGRA and MANKIND are registered trademarks and ENHANCED is a trademark of PMT Corporation. © Copyright 2011 PMT Corporation. All Rights Reserved

**PMT CORPORATION**  
QUALITY, INNOVATION AND CUSTOMER DRIVEN

**PMT CORPORATION**  
1500 Park Road  
Chanhassen, MN USA 55317  
www.pmtcorp.com

Call Toll Free ... 24 Hours a Day/365 Day a Year  
**1-800-MANKIND®** (626-5463)  
Phone 952-470-0866 » Fax 952-470-0865