

A FAMILY OF INNOVATION

For over three decades PMT® has been designing and producing quality electrodes and microwires for epilepsy procedures and research. Our unmatched signal clarity, product quality, and dependable customer service makes us an industry leader. With a large sales force spread throughout the U.S. and distributors around the world, we can be on-site to work with you to define products to match your specific requirements. We are continuously driving new product development and demanding the best quality to stay ahead of your needs. For more information or to schedule a trial, contact your local sales representative. For customers in the U.S., call **1-800-MANKIND**, for international customers call 1.952.470.0866.

FDA REGISTERED
AND APPROVED




PMT®
CORPORATION
QUALITY, INNOVATION AND CUSTOMER DRIVEN

1.800.MANKIND (626-5463) | WWW.PMTCORP.COM | SALES@PMTCORP.COM
FOR INTERNATIONAL CUSTOMERS | 1.952.470.0866 | INTLSALES@PMTCORP.COM
© 2011 PMT CORPORATION. ALL RIGHT RESERVED. ISO 13485: 2003 CERTIFIED

PRODUCT INFORMATION

PMT® CORTAC®


INTRACRANIAL ELECTRODES



PMT® CORTAC® ELECTRODES:

- » Unmatched signal clarity.
- » Reliable product quality.
- » Multiple design configurations.
- » Locally dedicated PMT product specialist.
- » Design your own configuration.
- » FDA registered and approved.

DESIGN CONFIGURATIONS



CONFIGURATIONS AVAILABLE:

- » Exits - grids are available with 1 to 8 exits depending on the number of contacts.
- » Contacts - available with 2.0mm, 3.0mm, or 4.5mm.
- » Material - contacts and wires available in platinum or stainless steel.
- » Thickness - strips and grids are available in .5mm, .76mm, or 1.01mm.
- » Surface - all grids and strips are available with enhanced or smooth surface.
- » Radiopaque - markers are available.
- » Wires - fully color coded wires for easy identification.

Local sales representatives are available 24 hours a day,
7 days a week, 365 days a year to assist you.
Contact PMT® Corporation for more information.

These products are protected by one or more of the following
U.S. patents: 6,162,101 & 5,961,909

PMT, CORTAC and MANKIND are registered trademarks of PMT Corporation.
© Copyright 2011 PMT Corporation. All Rights Reserved

FDA REGISTERED
AND APPROVED



PMT®
PMT CORPORATION
QUALITY, INNOVATION AND CUSTOMER DRIVEN

PMT CORPORATION
1500 Park Road
Chanhassen, MN USA 55317
www.pmtcorp.com

Call Toll Free ... 24 Hours a Day/365 Day a Year
1-800-MANKIND® (626-5463)
Phone 952-470-0866 » Fax 952-470-0865