

# CERVIMAX™ CTO

## *Cervical Thoracic Orthosis*

- CoolMax™ Lining
- All pediatric and adult sizes
- Effective patient immobilization
- MRI Compatible

*Provided to you exclusively by*

**PMT**®

*Quality, Innovation, Customer Driven*

1-800-MANKIND (626-5463)



Premium Quality

Maximum Benefits

Timeless Innovation

# Product Information

# CERVMAX™ CTO

## Description:

The CervMax™ CTO is designed to be used exclusively with the PMT® CervMax™ Cervical Collar.

Engineered for patients who require greater immobilization than a cervical collar alone can provide.

**Model 1625** CervMax™ CTO Brace with anterior only support bar. Includes extra set of replacement pads.

**Model 1625-0** CervMax™ CTO Brace with anterior and posterior support bar. Includes extra set of replacement pads.

**Model 1625-5** CervMax™ CTO Brace without extra set of replacement pads.

**Model 1625-10** Full set of replacement pads.

## Sizes

The CervMax™ CTO is available in adult and pediatric sizes to be used in conjunction with all models and sizes of the CervMax™ Cervical Collar.

## How Supplied

Each CervMax™ CTO Brace is provided standard with an extra set of replacement pads, patient care guide, and instruction manual.

## Special Features

- CoolMax™ fabric lines the inside of the CervMax™ CTO Brace.
- CoolMax™ dries faster than any other fabric. In moisture management and breathability tests vs. the competition, CoolMax™ is proven to have the fastest drying rate, quickest and most efficient movement of moisture away from the skin, and the best breathability.



## Our Commitment...Customer Satisfaction

PMT® is committed to providing you the highest quality products and the finest customer service.

Products specialists throughout the United States are available to serve you 24 hours a day, 7 days a week, 365 day a year. Please contact PMT® at 1-800-MANKIND (626-5463).

PMT and Mankind are registered trademarks and CervMax™ is a trademark of PMT® Corporation. Several Patents Pending. CoolMax™ is a DuPont trademark.

©2003 PMT® Corporation. All rights reserved.



PMT's CervMax™ CTO was engineered with your patients' comfort and stabilization in mind. Lightweight materials are MRI compatible, silicone free and latex free.