


More Fine Surgical  
Accessories from  
PMT.

SURGICAL SUCTION AND IRRIGATION  
**MACROVAC<sup>®</sup>, MICROVAC<sup>®</sup>  
& INVAC<sup>®</sup>**



*Benefits...*

**MacroVac<sup>®</sup>**

- Suction/irrigators - malleable tips, with various French sizes available

**MicroVac<sup>®</sup>**

- Ultra malleable stationary disposable suction accessory

**InVac<sup>®</sup>**

- Advanced valves for in-line suction control.
- Provides precision control simply by turning the outer collar to reduce or increase suction.

**PMT<sup>®</sup>**

Quality, Innovation, Customer Driven

# Product Information

## Description:

### **Model 3201 MacroVac - Designed for neurosurgery.**

- Malleable
- ESSC™ Enhanced Suction Control

Tip Styles: Straight, Teardrop, Medium, Ported

French Sizes: 3, 5, 7, 9, 12

Lengths: 3", 6", 9"

Finishes: Black Coat or Polished Stainless

### **\* Special Pricing available on neurological set of 12.**

When ordering, please indicate tip style, french size, length, and finish  
PMT is willing to customize any MacroVac to fit your needs

### **Model 3101 MicroVac - Ultramalleable stationary suction tips.**

French sizes: 3, 5, 7, 13, 22

22FR for laser smoke removal

When ordering, please indicate french size

### **Model 3102 InVac In-Line Suction Control**

**Model 3102-1** Surgical Size

**Model 3102-2** Microsurgical Size

**Model 3102-3** Macrosurgical Size

### **Model 3202-0 Ultrasoft™ Surgical Tubing for use with these products.**

Each tube is 6' in length. Available in:



Box (12 tubes)

Case (48 tubes)

## **Our Commitment...**

## **Customer Satisfaction**

PMT is committed to providing you the highest quality products and the finest customer service. Local sales representatives are available to serve you. Please contact PMT at 1-800-626-5463 or 952-470-0866.



PMT, MacroVac, MicroVac and InVac are registered trademarks and Ultrasoft and ESSC are trademarks of PMT Corporation.  
© Copyright 1998 PMT Corporation. All rights reserved.

**PMT®**  
*Quality, Innovation, Customer Driven*

PMT CORPORATION  
1500 Park Road  
Chanhassen, MN USA 55317

Call Toll Free  
**1-800-626-5463**  
Phone 952-470-0866 Fax 952-470-0865