



redefining beauty

Permark | UltraEssence®

The UltraEssence® provides precise digital needle speed which provides uniform pigment application throughout the procedure. The device can be set to manual mode by selecting one of four preset and programmable buttons. The foot pedal can be used to provide on demand speed variation. The UltraEssence is constructed of thick aluminum to create a solid, stable and secure device platform that will not slide or move during a procedure. The user screen has a touch panel for programming the machine to specific settings based on the procedure being performed and provides an LED screen for easily legible information including needle speed, procedure time and total machine time. The straight design precision hand piece accommodates all 53 shades of Permark pigments as well as all 12 needle variations. With so many features, options and variations the UltraEssence is the only micropigmentation device your spa or medical facility will ever require.

Permark®

1.800.MANKIND (626-5463) | WWW.PMTCORP.COM | SALES@PMTCORP.COM
FOR INTERNATIONAL CUSTOMERS | 1.952.470.0866 | INTLSALES@PMTCORP.COM

© 2010 PMT CORPORATION. ALL RIGHT RESERVED. PATENT PENDING. ISO 13485: 2003 CERTIFIED

Permark®

UltraEssence® includes the following:



UltraEssence™ includes Power Control Unit

- » Hospital grade medical casing
- » Programmable pre-set speed settings
- » Procedure time tracking
- » LED display screen
- » Powerful microprocessor can provide speeds up to 185 insertions per second Foot Pedal
- » Variable speed Precision Handpiece
- » Accomodates all 12 Permark® needles and 53 pigment shades
- » Straight design for maximum comfort and function
- » Reciprocating needle mechanism
- » Nose cone options for accurate depth control
- » Ability to sterilize by gas, autoclave or Sterrad®

PRODUCT INFORMATION

Permark, a division of PMT, offers the highest quality micropigmentation devices on the market complimented by the broadest selection of quality pigments, reusable and disposable needles. Our pigments are provided sterile in 53 shades for cosmetic and reconstructive procedures. Permark pigments are developed to be thick, creamy, concentrated formulas that provide no mess or splattering and require no water. They are formulated with purified FD&C ingredients and are hypoallergenic. The Permark needles are available in 12 variations for procedures ranging from fine lines such as eye liner to broad color application such as nipple areola reconstruction. Please contact your PMT/Permark product specialist for more information or call 1-800-PERMARK.

Sterile Pigment Color Options

P200050 Sandstone 1		P200162 Flesh 7		P200272 Brown/Black		P200442 Kajal		P201514 Autumn Red	
P200060 Sandstone 2		P200172 Flesh 8		P200222 Black		P200422 Taupe		P201515 Desert Rose	
P200070 Sandstone 3		P200182 Flesh 9		P200232 Dark Gray		P201501 Red		P201517 Night Red	
P200090 Flesh .5		P200192 Flesh 10		P200301 Navy Blue		P201500 Mauve		P201520 Chocolate	
P200100 Flesh 1		P200202 Flesh 11		P200303 Slate Blue		P201502 Berry Red		P201521 Dark Chocolate	
P200101 Flesh 1.5		P200212 Flesh 12		P200238 Pewter		P201503 Mocha		P200485 White Toner	
P200102 Flesh 2		P200432 Topaz		P200294 Charcoal Green		P201507 Brick		P200486 Brown Toner	
P200112 Flesh 2.5		P200282 Light Brown		P200296 Green		P201505 Coral		P200487 Yellow Toner	
P200132 Flesh 4		P200252 Med. Brown		P200234 Medium Gray		P201506 Terra Cotta		P200488 Red Toner	
P200142 Flesh 5		P200292 Umber		P200452 Blonde		P201513 Plum			
P200152 Flesh 6		P200242 Dark Brown		P200412 Blond K21		P201510 Blush			

