

PMT sEEG ELECTRODES

2102 SERIES



PMT sEEG's are less invasive due to a smaller diameter (0.8mm).

sEEG's are offered with removable stylet or permanent stylet.

PMT sEEG's are available as nonmagnetic.

Our anchor bolts attach to the anchor bolt driver and defies gravity during insertion. They are also designed to reduce rotation of the sEEG electrode during tightening of the cap.

PMT's sEEG's are stocked on the shelf ready to ship.

PMT's sEEG procedure is shorter due to minimized steps for anchor bolt and electrode placement.

PMT's instruments all fit through a 2.4mm OD bushing allowing the elimination of reducing tubes.

PMT offers a blank cap to stop CFS leakage through the anchor bolts before sEEG implantation.

PMT offers stereotactic placement kits to work with nearly every major stereotactic frame/system.

PMT sEEG Electrodes are available with 8, 10, 12, 14 or 16 contacts.




Image provided by Memorial Hermann and George Craig Photography

PRODUCT INFORMATION

PMT sEEG ELECTRODES

- » PMT sEEG Electrodes are available with 8, 10, 12, 14 or 16 contacts.
- » All models are available with either a permanent stylet or a removable stylet
- » All models are 0.8mm in diameter and work with any stereotactic system

MODEL NUMBER	DESCRIPTION	STYLET	CABLE NUMBER
2102-08-091	8 contacts, 26.5mm Recording depth	Removable	2120-08-01B
2102-08-094	8 contacts, 26.5mm Recording depth	Permanent	2120-08-01B
2102-10-091	10 contacts, 32.5mm Recording depth	Removable	2120-10-01B
2102-10-094	10 contacts, 32.5mm Recording depth	Permanent	2120-10-01B
2102-12-091	12 contacts, 40.5mm Recording depth	Removable	2120-12-01B
2102-12-094	12 contacts, 40.5mm Recording depth	Permanent	2120-12-01B
2102-14-091	14 contacts, 47.5mm Recording depth	Removable	2120-14-01B
2102-14-094	14 contacts, 47.5mm Recording depth	Permanent	2120-14-01B
2102-16-091	16 contacts, 54.5mm Recording depth	Removable	2120-16-01B
2102-16-092	16 contacts, 61.5mm Recording depth	Removable	2120-16-01B
2102-16-093	16 contacts, 68.5mm Recording depth	Removable	2120-16-01B
2102-16-094	16 contacts, 54.5mm Recording depth	Permanent	2120-16-01B
2102-16-095	16 contacts, 61.5mm Recording depth	Permanent	2120-16-01B
2102-16-096	16 contacts, 68.5mm Recording depth	Permanent	2120-16-01B
2102-16-097	16 contacts, 80mm Recording depth	Permanent	2120-16-01B
2102-16-098	16 contacts, 90mm Recording depth	Permanent	2120-16-01B

All sEEG placement kits shown are for use with a 2.4mm drill bit.

STEALTH sEEG KIT

MODEL #2103-41

BRW/CRW sEEG KIT

MODEL #2103-35

ROSA sEEG KIT

MODEL #2103-32

VERTEK sEEG KIT

MODEL #2103-43

LEKSELL sEEG KIT

MODEL #2103-37

BRAINLAB sEEG KIT

MODEL #2103-109

© Copyright 2014 PMT® Corporation. All Rights Reserved